

WIDS - AI For Healthcare

WEEK 03:

TASK 01:


Dataset: <https://universe.roboflow.com/cxrdataset/skin-ga5ww>

- **Type of images:** Skin Images
- **Number of samples:** 7206 Images
- **Classes:** This dataset has 4 classes.
 - I. **Basal Cell Carcinoma:** a type of skin cancer.
 - II. **Melanoma:** a serious form of skin cancer.
 - III. **Nevus:** usually a benign mole.
 - IV. Unlabeled
- **Dataset Splits:** The dataset is divided into training, validation, and test sets.
 - I. Training set: 5044 images
 - II. Validation set: 1442 images
 - III. Test set: 720 images
- **Imbalance :**
 - I. Some classes have more images than others.
 - II. Skin tone variation.
 - III. Lighting differences.


TASK 02:

- **Model Used:**
 - **Model:** YOLOv8 Classification Model
 - **Framework:** Ultralytics
 - **Training Platform:** Google Colab
 - **Hardware Used:** GPU (when enabled)YOLOv8 was chosen because it provides fast training and good performance for image classification tasks.
- **Training Results:**
 - **Number of Epochs:** 20
 - **Image Size:** 224 × 224
 - **Training / Validation Accuracy:** 98.3%
 - **Training Loss:** 0.0639


train/loss




metrics/accuracy_top1





val/loss



metrics/accuracy_top5



- Confusion Matrix** : The confusion matrix shows that the model correctly classifies most images.
 Minor confusion exists due to visual similarity between skin conditions.
-



- **Sample Predictions :** Sample predictions show correct classification of Basal Cell Carcinoma and Melanoma images.
- **Brief Observations :**
 - The model achieves high accuracy on available classes.
 - Limited class availability affected multi-class training.
 - YOLOv8 provides fast and efficient classification.
 - Dataset quality strongly impacts performance.
 - GPU acceleration significantly reduces training time.